

A GloMax[®]-Multi Jr Method for Coomassie Plus— The Better Bradford[™] Assay Kit

INTRODUCTION

The GloMax[®]-Multi Jr in combination with the Pierce Coomassie Plus Assay Kit provides a convenient procedure for quantifying protein. When the Coomassie dye binds protein, the absorption maximum shifts from 465 nm to 595 nm. The Absorbance Module detects as little as 25 µg of BSA in 1.5 mL of Coomassie Plus Reagent.




Figure 1. Bradford Assay was performed on the GloMax[®]-Multi Jr using the Absorbance Module, the 600 nm filter, and the Coomassie Plus reagent. Each point represents the average of replicate samples after background subtraction ($n=2$).

MATERIALS

- GloMax[®]-Multi Jr Multimode Reader
- Absorbance Module
- 600 nm filter paddle
- 10 x 10 mm disposable methacrylate cuvettes
- Coomassie Plus—The Better Bradford Assay Kit (Pierce Cat.# 23225, 23227)

PREPARATION

Note: Store assay reagent at 4°C. Unopened vials of bovine serum albumin (BSA) standard may be stored at room temperature.

1. BSA Standard Curve

Prepare a serial dilution of BSA that covers the range for your samples. For example, create a twofold dilution series from 1000 µg/mL to 125 µg/mL. Make sure to include a blank solution (diluent only) in your standard curve preparation.

Note: Prepare the dilution series in the same diluent as the samples for best results.

2. Coomassie Plus Reagent

- Determine the total volume of reagent required. Each sample and standard requires 1.5 mL of reagent.
- Gently invert the bottle of Coomassie Plus Reagent solution to mix the solution before removing the amount necessary for the assay.
- Equilibrate the reagent to room temperature before use.
- **Note:** Coomassie Plus Reagent may form dye-dye and dye-protein aggregates when left undisturbed. Gentle mixing dissolves the aggregates. For best results, mix the reagent before dispensing and again before measuring absorbance.

3. Samples

- For each standard or sample, pipette 50 µL into an individual 10 x 10 mm methacrylate cuvette.
- Add 1.5 mL of the Coomassie Plus reagent to each cuvette.
- Incubate the samples and standards for 10 minutes at room temperature.

4. Instrument Setup

- Power OFF the GloMax[®]-Multi Jr. Install the Absorbance Module into the sample compartment according to *Technical Manual*.
- Insert the 600 nm filter paddle in the Absorbance Module.
- Turn ON the GloMax[®]-Multi Jr, and use the touch screen to choose the Absorbance operation mode.
- Touch “Calibrate,” and use the black cuvette to set the GloMax[®]-Multi Jr to calibrate the zero (dark) reading.
- Use a cuvette containing 2 mL of appropriate buffer to calibrate the baseline (100% transmittance) reading.
- Touch “OK” to accept the calibrations and return to the “Home” screen.

SAMPLE ANALYSIS

- Insert the sample or standard into the Absorbance Module, and touch “Measure Absorbance” to begin measurement.
- Record the results in Absorbance units (Ab).
- Use a standard curve to determine the protein concentration of each unknown sample. A four-parameter (quadratic) or best-fit curve provides the best accuracy.

GloMax is a trademark of Promega Corporation.
Coomassie-Plus – The Better Bradford is a registered trademark of Pierce Biotechnology, Inc.

CONTACT INFORMATION

Toll-Free: (800) 356-9526
Fax: (800) 356-1970

www.promega.com

Email: custserv@promega.com

Mailing Address:

Promega Corporation
2800 Woods Hollow Rd.
Madison, WI 53711 USA